STERILE EO

CAUTION: Federal (U.S.A.) Law restricts this device to sale by or on the order of a physician.

Pressure Products

Customer Service
1861 N. Gaffey Street, Suite B
San Pedro, CA 90731 - USA
Tel +1-310-547-4973
Fax +1-310-547-4760

www.pressure-products.com

Merit Medical Systems, Inc.
65 Great Valley Parkway
Malvern, PA 19355 - USA
Tel +1-801-253-1600
USA Customer Service:
Tel +1-800-356-3748

Merit Medical Ireland, Ltd
Parkmore Business Park West
Galway, Ireland
European Customer Service:
Tel +31 43 358 8222

0344

For U.S.-California Only.

Proposition 65: A state of California voter initiative, requires the following notice:



WARNING: This product and its packaging have been sterilized with ethylene oxide. This packaging may expose you to ethylene oxide, a chemical known to the state of California to cause cancer or birth defects or other reproductive harm.

The H. Lee SafeSheath[®] hemostatic splittable sheath is protected by Patent Insurance. USA Patent 7,462,184. European Patent EP1360972. Other USA and Worldwide Patents pending. SafeSheath and CSG are registered trademarks of Pressure Products Medical Supplies, Inc.

ACCESSING THE HEART OF INNOVATION

Technical Manual

en-Directions for Use

This device is intended for one time use only. Read instructions prior to use.

Indications

For the introduction of various types of pacing or defibrillator leads and catheters.

Contraindications

Use of the Coronary Sinus Guide / Lateral Vein Introducer (CSG/LVI) systems is contraindicated for the following:

- Patients with an existing or possible occlusion of the coronary vessels or unsuitable anatomy of the coronary veins
- Patients with active systemic infection

Possible Negative Side Effects / Adverse Events

Coronary Sinus Guide / Lateral Vein Introducer (CSG, LVI) systems should be used by physicians familiar with percutaneous catheter introduction. Complications which may be associated with the use of catheter introducer systems include, but are not limited to, the following:

- Air embolism
- Allergic reaction to contrast media
- Arterial wall damage
- Bleeding
- Cardiac arrhythmias
- Cardiac tamponade
- Chronic nerve damage
- Damage to the heart valves
- Hematoma at the puncture site
- Infection
- Local tissue response, fibrotic tissue formation
- Myocardial damage
- Myocardial infarction
- Plaque dislodgement
- Pneumothorax
- Stroke and death
- Thrombus formation/emboli
- Vascular occlusion
- Vascular spasm
- Venous or cardiac perforation

1

Warnings

- This product is sensitive to light. Do not use if stored outside the protective outer carton. Store in a cool, dark, and dry place.
- Infusion through the side port can be done only after all air is removed from the unit. Improper use of the transvalvular insertion tool (TVI) can cause air embolism and back bleeding.
- Do not use this device in patients who cannot be appropriately anticoagulated. When tested in non-anticoagulated sheep, this device has shown thrombus formation; however, heparinized studies alleviated the concern.

Precautions

- Do not alter this device in any way.
- **Single Use Devices:** This single-use product is not designed or validated to be reused. Reuse may cause a risk of cross-contamination, affect the measurement accuracy, system performance, or cause a malfunction as a result of the product being physically damaged due to cleaning, disinfection, re-sterilization, or reuse.
- Federal law (U.S.A.) restricts this device to sale by or on the order of a physician.
- Aspiration and saline flushing of the sheath, dilator, and valve should be performed to help minimize the potential for air embolism and clot formation.
- Indwelling introducer sheaths should be internally supported by a catheter, pacing lead, or dilator.
- Dilators, catheters, and pacing leads should be removed slowly from the sheath. Rapid removal may damage the valve members resulting in blood flow through the valve. Never advance or withdraw guide wire or sheath when resistance is met. Determine cause by fluoroscopy and take remedial action.
- When injecting or aspirating through the sheath, use the side port only.
- When using the transvalvular insertion tool (TVI), lead size may not exceed 6.2F.
- When using the TVI always keep the exposed proximal end covered to prevent air embolization and back bleeding.

en Directions for Use	1
da Brugsvejledning	2
de Gebrauchsanweisung	3
el Οδηγίες χρήσης	4
es Instrucciones de uso	5
fi Käyttöohje	6
fr Mode d'emploi	6
it Istruzioni per l'uso	7
ja 使用説明書	8
nl Gebruiksaanwijzing	9
no Bruksanvisning	10
pt Instruções de Utilização	11
sv Bruksanvisning	12

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