








Braided Sheath Cutter

For use with SafeSheath® CSG Braided Introducers

Distributed by

Pressure Products, Inc.

Customer Service
1861 N. Gaffey Street Suite B
San Pedro, CA 90731 USA
1.310.547.4973 1.800.600.4973
Fax 1.310.547.4760
www.pressure-products.com

Manufactured by

Thomas Medical Products, Inc.

65 Great Valley Parkway
Malvern, PA 19355 USA
1.610.296.3000 Fax 1.610.296.4591

Authorized European Representative

Vital Signs, Ltd.

13-14 Eldon Way
Lineside Industrial Estate
Littlehampton, West Sussex, UK

The SafeSheath®, hemostatic splittable sheath is protected by Patent Insurance. USA Patents 5125904 and 5409463. Canadian Patent 2,112,394. Worldwide Patents pending. The Braided Sheath Cutter is covered by USA Patent 6,497,681. Other patents pending.

SafeSheath and CSG are a registered trademarks of Pressure Products, Inc.

Directions for Use

To aid in removal of a braided sliceable introducer after it has been used to place a pacemaker lead or catheter.

All Cutter materials are bio-compatible, latex-free and PVC-free.

Precautions

Do not alter this device in any way.

Federal law (U.S.A.) restricts this device to sale by or on the order of a physician.

As specified in the following chart, use an introducer sheath and pacing lead within the parameters listed in the chart below.

Largest Introducer Size	Largest Lead Size
4.0mm/12Fr	2.0mm/6Fr

Suggested Directions for Use

The Cutter is to be used after pacing lead or catheter is in the desired location and the sheath introducer is ready to be removed. It is suggested a lead stylet or stiffening wire be retained in the pacing lead during sheath removal.

■ If cutting a SafeSheath CSG Braided Series Sheath:

1. Open package and place contents on sterile field.
2. It is suggested that the braided sheath be withdrawn to the furthest extent possible prior to splitting the sheath valve and cutting the sheath.

NOTE: The lead tip and lie of the lead along the RA floor should be carefully observed fluoroscopically as the braided sheath is withdrawn and cut, to avoid inadvertent lead dislodgment.

3. Manually split the handles of the sheath valve. One half of the valve will remain attached to the braided sheath. Dispose of the other valve section. (see Figure 1)

4. Position the distal 'C' section of the Cutter under the sheath tubing, distal to the sheath hub. (see Figure 2)

5. Place the blade hook into the 'V' notch of the tubing at the proximal end of the sheath. (see Figure 3)

6. Thread the pacing lead over the notch at the back end of the Cutter and compress the lead with your thumb, while holding the Cutter in a 'pistol grip' fashion. (see Figure 4)

7. While holding the Cutter and the lead with one hand, slowly withdraw the braided sheath over the Cutter blade with the other hand. As the sheath is withdrawn it is cut away from the lead in the process. (see Figure 4)

■ If cutting a braided Lateral Vein Introducer (LVI) within a CSG sheath:

1. Open package and place contents on sterile field.
2. With a stabilizing stylet or stiffening wire in place, and while holding the lead in place, under fluoroscopic control, withdraw the LVI into the distal end of the SafeSheath CSG. (see Figure 5)

NOTE: The lead tip should be carefully observed fluoroscopically as the LVI sheath is withdrawn and cut, to avoid inadvertent lead dislodgment. The outer CSG sheath should be stabilized by an assistant grasping the CSG distal to the CSG hub as the LVI is withdrawn and cut to prevent inadvertent dislodgment of the CSG outer sheath from the mid CS.

3. Manually split the handles of the LVI sheath valve. One half of the valve will remain attached to the braided LVI. Dispose of the other valve section. (see Figure 6)

4. Position the distal 'C' section of the Cutter under the LVI sheath tubing distal to the sheath hub. (see Figure 7)

5. Place the blade hook into the 'V' notch of the LVI tubing at the proximal end of the LVI guide. (see Figure 3)

6. While maintaining the LVI Guide, pacing lead and CSG Sheath in position, hold the Cutter in 'pistol grip' fashion and cut the LVI by advancing the Cutter toward the CSG sheath until the distal 'C' portion of the Cutter engages the CSG hub. (see Figure 8)

NOTE: Make sure the pacing lead, LVI guide and CSG sheath are not advanced or withdrawn with the Cutter during this step of the procedure.

7. After the Cutter has engaged the CSG hub, thread the pacing lead over the notch at the back end of the Cutter and compress the lead with your thumb. (see Figure 9)

8. While holding the Cutter engaged with the CSG hub and lead with one hand, slowly withdraw the LVI sheath over the Cutter blade with the other hand. As the LVI is withdrawn it is cut away from the lead in the process. (see Figure 9)

NOTE: The outer CSG sheath should continue to be stabilized distal to the CSG hub as the LVI is withdrawn and cut to prevent inadvertent dislodgment of the CSG outer sheath from the mid CS.

9. After the LVI is cut and removed the CSG is removed as described above under the heading "If cutting a SafeSheath CSG Braided Series Sheath," or in the appropriate DFU.

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Føderale lover (U.S.A.) begrenser salg av denne anordning til eller på anbefaling av en lege.