
CPS®

H. Lee
SafeSheath®**SEALING ADAPTER****STERILE EO****For U.S.-California Only.****Proposition 65**, a State of California voter initiative, requires the following notice:**WARNING:** This product and its packaging have been sterilized with ethylene oxide. This packaging may expose you to ethylene oxide, a chemical known to the state of California to cause cancer or birth defects or other reproductive harm.**CAUTION: Federal (U.S.A.) Law restricts this device to sale by or on the order of a physician.****Cardiac Positioning System (CPS®)**

en	Cardiac Positioning System (CPS®) Sealing Adapter
cs	<i>Systém pro zavádění srdečních elektrod (CPS®) Těsnicí Adaptér</i>
da	<i>Cardiac Positioning System (CPS®) Forseglings adapter</i>
de	<i>Cardiac Positioning System (CPS®) Abdichtungsadapter</i>
el	<i>Cardiac Positioning System (CPS®) Προσαρμογέας Σφράγιση</i>
es	<i>Sistema de Posicionamiento Cardíaco (CPS®) Adaptador de Sellado.</i>
et	<i>CPS® (Südamepositsioneerimissüsteem) Tihendusadapter</i>
fi	<i>CPS®-järjestelmä (Cardiac Positioning System) Sulkuliitin</i>
fr	<i>Cardiac Positioning System (CPS®) Raccord d'étanchéité</i>
hu	<i>CPS® (Kardiális pozícionáló rendszer Tömítőadapter</i>
it	<i>Sistema di Posizionamento Cardíaco (CPS®) Adattatore a tenuta</i>
ja	<i>Cardiac Positioning System (CPS®) 密閉アダプタ</i>
ko	<i>CPS® (Cardiac Positioning System) 실링 어댑터</i>
lt	<i>CPS® Širdies padėties nustatymo sistema Sandarinimo adapteris</i>
nl	<i>Cardiac Positioning System (CPS®) Sealingadapter</i>
pl	<i>System pozycjonowania serca (CPS®) Nasadka uszczelniająca</i>
pt	<i>Sistema de Posicionamento Cardíaco (CPS®) Adaptador de vedação</i>
ru	<i>Система позиционирования сердца (CPS®) Герметичный переходник</i>
sk	<i>Systém pre umiestňovanie intrakardiálnych elektrod (CPS®) Tesniaci adaptér</i>
sv	<i>CPS® (Cardiac Positioning System) Tätningsadapter</i>
zh	<i>心脏定位系统 (CPS®) 密封接头</i>

EN	Directions for Use	2
CS	Návod k použití	2
DA	Brugsanvisning	3
DE	Gebrauchsanleitung	3
EL	Οδηγίες χρήσης	3
ES	Instrucciones para el uso	3
ET	Kasutusjuhised	4
FI	Käyttöohjeet	4
FR	Mode d'emploi	4
HU	Használati utasítás	4
IT	Istruzioni per l'uso	5
JA	使用説明	5
KO	사용 방법	5
LT	Naudojimo instrukcijos	5
NL	Gebruiksaanwijzing	6
PL	Instrukcja obsługi	6
PT	Instruções de Utilização	6
RU	Инструкции по использованию	6
SK	Návod na použitie	7
SV	Bruksanvisning	7
ZH	使用说明	7

St. Jude Medical – Cardiac Rhythm Management Division
15900 Valley View Court
Sylmar, CA 91342 USA
Tel +1-818.362.6822
Tel 1.800.722.3774 (US only)

Pressure Products
Customer Service
1861 N. Gaffey Street, Suite B
San Pedro, CA 90731 - USA
Tel +1-310-547-4973
Fax +1-310-547-4760
www.pressure-products.com

CE
0344

Merit Medical Systems, Inc.
65 Great Valley Parkway
Malvern, PA 19355 - USA
Tel +1-801-253-1600
USA Customer Service:
Tel +1-800-356-3748

EC REP








Merit Medical Ireland, Ltd
Parkmore Business Park West
Galway, Ireland
European Customer Service:
Tel +31 43 358 82 22

The H. Lee SafeSheath® hemostatic splitable sheath is protected by Patent Insurance. USA Patent 8,048,061. Canadian Patent 2,493,189. European Patent EP1545677. Japanese Patent 4705779. Other Worldwide Patents pending. SafeSheath is a registered trademark of Pressure Products Medical Supplies, Inc. CPS® is a registered trademark of St. Jude Medical, Inc.

ST. JUDE MEDICAL®
MORE CONTROL. LESS RISK.

Pressure Products
ACCESSING THE HEART OF INNOVATION

Technical
Manual**(EN)-Directions for Use**

This device is intended for one time use only. Read instructions prior to use.

INDICATIONS

For the introduction of various types of pacing leads and catheters.

WARNINGS

Infusion through the sideport can be done only after all air is removed from the unit.

PRECAUTIONS

- Do not alter this device in any way.
- **Single Use Devices:** This single-use product is not designed or validated to be reused. Reuse may cause a risk of cross-contamination, affect the measurement accuracy, system performance, or cause a malfunction as a result of the product being physically damaged due to cleaning, disinfection, re-sterilization, or reuse.
- Aspiration and saline flushing of the CPS® / SafeSheath® Sealing Adapter should be performed to help minimize the potential for air embolism and clot formation.
- Dilators, catheters, pacing leads and other leads should be removed slowly from the CPS® / SafeSheath® Sealing Adapter. Rapid removal may damage the valve members resulting in blood flow through the valve.
- When injecting or aspirating through the CPS® / SafeSheath® Sealing Adapter, use the sideport only.

**USE STERILE TECHNIQUE
A SUGGESTED PROCEDURE**

1. Align the base of the hemostasis valve assembly over the proximal end of the sheath hub so the sideport tubing is perpendicular to the T-handle. (see figure 1)
2. Firmly press the base of the hemostasis valve into the proximal hub of the sheath in the middle of the T-handle Hub by aligning the slots on the Adapter with the locking tabs on the CPS® Hub (see figure 2) and push the Adapter into the CPS® Hub until the top of the Sealing Adapter elastic plug is flush with top of the Hub, taking care not to split the sheath hub or T-handle.
3. Rotate the hemostasis valve on the sheath hub so the sideport tubing is positioned over either side of the T-handle, being sure the two raised flanges on the sides of the hemostasis valve engage the locking tabs at the base of each side of the T-handle. (see figure 3)
4. Flush all components of the Catheter Delivery System with heparinized saline solution.
5. When removing the outer catheter, disconnect the hemostasis valve from the outer guide catheter hub by breaking the outer guide catheter's T-handle and tearing it away from the hemostasis valve. (see figure 4 or 5)
6. Sharply snap the tabs of the CPS® / SafeSheath® Sealing Adapter down and away to separate the valve from the lead. (see figure 6 or 7)

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