**STERILE | EO**

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


## en-Instructions for Use

This device is intended for one time use only. Read instructions prior to use.

### Indications

For the introduction of various types of pacing leads and catheters.

### Warnings

- This product is sensitive to light. Do not use if stored outside the protective outer carton. Store in a cool, dark, and dry place.
- Infusion through the side port can be done only after all air is removed from the unit.

### Precautions

- Do not alter this device in any way.
- Single Use Devices:** This single-use product is not designed or validated to be reused. Reuse may cause a risk of cross-contamination, affect the measurement accuracy, system performance, or cause a malfunction as a result of the product being physically damaged due to cleaning, disinfection, re-sterilization, or reuse.
- Aspiration and saline flushing of the sheath, dilator, and valve should be performed to help minimize the potential for air embolism and clot formation.
- Indwelling introducer sheaths should be internally supported by a catheter, electrode, or dilator.
- Dilators, catheters, and pacing leads should be removed slowly from the sheath. Rapid removal may damage the valve members resulting in blood flow through the valve. Never advance or withdraw guide wire or sheath when resistance is met. Determine cause by fluoroscopy and take remedial action.
- When injecting or aspirating through the sheath, use the side port only.

### Use Sterile Technique A suggested procedure:

- Peel open package and place contents on sterile field.
- Prep skin and drape in area of anticipated venipuncture as desired.
- Distend the subclavian vein or internal jugular vein. The subclavian vein is difficult to locate unless it is distended by raising the patient's legs to a 45° angle or by using the Trendelenburg position. The vein will be much easier to locate if the patient is well hydrated.
- Insert needle into vessel. The needle position should be verified by observing venous blood return. (Puncture the vein, as medially as possible, by entering just lateral to the ligament which joins the clavicle and the first rib.)
- The angle of the needle should be adjusted depending on the patient's build: shallow in a thin person, deeper in a heavyset person. Use an 18 gauge needle, 7cm (2-3/4 in.) long.
- Aspirate the puncture needle using the 10cc syringe.
- Remove the syringe and insert soft tip of guide wire through the introducer needle into the vessel. Advance guide wire to required depth. Leave an appropriate amount of guide wire exposed.

**At no time should the guide wire be advanced or withdrawn when resistance is met.** Determine the cause of resistance before proceeding. Fluoroscopic verification of the guide wire's entrance into the superior vena cava and right atrium is suggested.

- Guide wire in place and remove introducer needle. Do not withdraw the guide wire back into the cannula as this may result in separation of the guide wire. The cannula should be removed first.
- Insert vessel dilator into sheath until the dilator cap

folds over valve housing and secures the dilator onto sheath assembly.

- Thread the dilator/sheath assembly over the guide wire.
- Advance the dilator and sheath together with a twisting motion over the guide wire and into the vessel. Fluoroscopic observation may be advisable. Attaching a clamp or hemostat to the proximal end of the guide wire will prevent inadvertently advancing the guide wire entirely into the patient.
- Once assembly is fully introduced into the venous system, separate the dilator cap from the sheath valve housing by rocking the dilator cap off the hub. (see figure 1)
- Slowly retract the guide wire and dilator, leaving the sheath in position. The hemostasis valve will reduce the loss of blood and the inadvertent aspiration of air through the sheath.
- Aspirate all air from the sheath valve assembly by using a syringe connected to the side port. Flush the introducer through the side port. If the introducer is to remain in place during lead positioning and testing, flush the introducer via the side port periodically with saline is advised.
- Introduce pacemaker lead or catheter through the hemostasis valve/sheath and advance it into position.
- Flush sheath with 5cc of saline immediately before peeling sheath away in order to minimize backbleeding.
- Withdraw sheath and valve over the lead or catheter and from the vessel, while keeping the lead in place.
- Sharply snap the tabs of valve housing in a plane perpendicular to the long axis of the sheath to split the valve and peel sheath apart while withdrawing from the vessel. (see figure 2)
- A retained guide wire technique may be used for dual lead implantation. The valve will remain hemostatic with both the lead and guide wire inserted through it; however, extreme caution must be used when manipulating the lead in order to prevent accidental advancement of the guide wire into the patient. Again, a hemostat attached to the guide wire proximal end is advised.

## TVI Directions for Use

This device is intended for one time use only. Read instructions prior to use.

**NOTE:** A transvalvular insertion (TVI) tool is to be used at the physician's discretion to open the valve of a tear-away hemostatic introducer for ease of lead placement.

### Indications

For the introduction of various types of pacing or defibrillator leads and catheters.

### Warnings

Improper use of a Transvalvular Insertion Tool (TVI) can cause air embolism and back bleeding.

### Precautions

When the TVI is inserted into the tear-away hemostatic introducer valve housing all hemostasis is lost and the risk of air embolization and back bleeding exists. Keep the proximal exposed end of the TVI covered with your thumb while the TVI is in use.

### Use of a Transvalvular Insertion Tool (TVI) with a tear-away hemostatic introducer

*A suggested procedure:*

- Hold your thumb over the proximal exposed opening of the TVI to minimize the potential for air embolization or back bleeding.
- Insert the distal end of the TVI into the tear-away hemostatic introducer valve by gently pushing the TVI through the valve just until the distal end of the TVI contacts the interior of the introducer hub. (see figure 3)

**Ensure the TVI has passed entirely through the valve and the valve is fully dilated before attempting to introduce the lead through the TVI. Do not over-insert the TVI, if the TVI is over-inserted it can cause damage to the tool and impact lead passage.**

- Advance the lead through the TVI and into the sheath. If the lead is difficult to advance, stop and remove the lead and TVI, and replace with a new TVI.
- As soon as the lead body or shocking coils are resting inside the sheath and past the TVI, pull the TVI back out of the sheath valve housing.
- Once the TVI has been withdrawn from the introducer, aspirate the introducer through the sideport until any air which may have entered the sheath during the procedure is removed, and again flush with heparinized saline.
- The TVI may be separated and removed, or it may be temporarily left resting on the body of the lead to assist in further repositioning. Once the lead is in position, the TVI must be removed from the lead.

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SafeSheath is a registered trademark of Pressure Products Medical Supplies, Inc.