**STERILE EO**

**For U.S.-California Only.**  
**Proposition 65:** A state of California voter initiative, requires the following notice:

**WARNING:** This product and its packaging have been sterilized with ethylene oxide. This packaging may expose you to ethylene oxide, a chemical known to the state of California to cause cancer or birth defects or other reproductive harm.

#### Pressure Products

Customer Service  
 1861 N. Gaffey Street, Suite B  
 San Pedro, CA 90731 - USA  
 Tel +1-310-547-4973  
 Fax +1-310-547-4760  
[www.pressure-products.com](http://www.pressure-products.com)

**Merit Medical Systems, Inc.**  
 65 Great Valley Parkway  
 Malvern, PA 19355 - USA  
 Tel +1-801-253-1600  
 USA Customer Service:  
 Tel +1-800-356-3748

**EC REP**

**Merit Medical Ireland, Ltd**  
 Parkmore Business Park West  
 Galway, Ireland  
 European Customer Service  
 Tel +31 43 358 82 22

**CE**

**0344**

**Pressure®  
Products**  
ACCESSING THE HEART OF INNOVATION

**CAUTION: Federal (U.S.A.) Law restricts this device to sale by or on the order of a physician.**

## Technical Manual

### en-Directions for Use

This device is intended for one time use only. Read instructions prior to use.

#### Indications

For the introduction of various types of pacing leads and catheters.

#### Warnings

- This product is sensitive to light. Do not use if stored outside the protective outer carton. Store in a cool, dark, and dry place.
- Infusion through the side port can be done only after all air is removed from the unit.

#### Precautions

- Do not alter this device in any way.
- **Single Use Devices:** This single-use product is not designed or validated to be reused. Reuse may cause a risk of cross-contamination, affect the measurement accuracy, system performance, or cause a malfunction as a result of the product being physically damaged due to cleaning, disinfection, re-sterilization, or reuse.
- Aspiration and saline flushing of the sheath, dilator, and valve should be performed to help minimize the potential for air embolism and clot formation.
- Indwelling introducer sheaths should be internally supported by a catheter, electrode, or dilator.
- Dilators, catheters, and pacing leads should be removed slowly from the sheath. Rapid removal may damage the valve members resulting in blood flow through the valve. Never advance or withdraw guide wire or sheath when resistance is met. Determine cause by fluoroscopy and take remedial action.
- When injecting or aspirating through the sheath, use the side port only.


#### Use Sterile Technique A suggested procedure:

1. Peel open package and place contents on sterile field.
2. Prep skin and drape in area of anticipated venipuncture as desired.
3. Distend the subclavian vein or internal jugular vein. The subclavian vein is difficult to locate unless it is distended by raising the patient's legs to a 45° angle or by using the Trendelenburg position. The vein will be much easier to locate if the patient is well hydrated.
4. Insert needle into vessel. The needle position should be verified by observing venous blood return. (Puncture the vein, as medially as possible, by entering just lateral to the ligament which joins the clavicle and the first rib.)
5. The angle of the needle should be adjusted depending on the patient's build: shallow in a thin person, deeper in a heavyset person. Use an 18-gauge needle, 7 cm (2-3/4 in.) long.
6. Aspirate the puncture needle using the 12cc syringe.
7. Remove the syringe and insert soft tip of guide wire through the introducer needle into the vessel. Advance guide wire to required depth. Leave an appropriate amount of guide wire exposed.


**At no time should the guide wire be advanced or withdrawn when resistance is met.** Determine the cause of resistance before proceeding. Fluoroscopic verification of the guide wire's entrance into the superior vena cava and right atrium is suggested.

8. Hold guide wire in place and remove introducer needle. Do not withdraw the guide wire back into the cannula as this may result in separation of the guide wire. The cannula should be removed first.
9. Insert vessel dilator into sheath until the dilator cap folds over valve housing and secures the dilator onto sheath assembly.
10. Thread the dilator/sheath assembly over the guide wire.
11. Advance the dilator and sheath together with a twisting motion over the guide wire and into the vessel. Fluoroscopic observation may be advisable. Attaching a clamp or hemostat to the proximal end of the guide wire will prevent inadvertently advancing the guide wire entirely into the patient.
12. Once assembly is fully introduced into the venous system, separate the dilator cap from the sheath valve housing by rocking the dilator cap off the hub. (see figure 1)
13. Slowly retract the guide wire and dilator, leaving the sheath in position. The hemostasis valve will reduce the loss of blood and the inadvertent aspiration of air through the sheath.
14. Aspirate all air from the sheath valve assembly by using a syringe connected to the side port. Flush the introducer through the side port. If the introducer is to remain in place during lead positioning and testing, flushing the introducer via the side port periodically with saline is advised.
15. Introduce pacemaker lead or catheter through the hemostasis valve/sheath and advance it into position.
16. Flush sheath with 5cc of saline immediately before peeling sheath away in order to minimize backbleeding.
17. Withdraw sheath and valve over the lead or catheter and from the vessel, while keeping the lead in place.
18. Sharply snap the tabs of valve housing in a plane perpendicular to the long axis of the sheath to split the valve and peel sheath apart while withdrawing from the vessel. (see figure 2)
19. A retained guide wire technique may be used for dual lead implantation. The valve will remain hemostatic with both the lead and guide wire inserted through it; however, extreme caution must be used when manipulating the lead in order to prevent accidental advancement of the guide wire into the patient. Again, a hemostat attached to the guide wire proximal end is advised.

**Figure 1**

**Figure 2**

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The H. Lee SafeSheath® hemostatic splittable sheath is protected by Patent Insurance. USA Patent 5,409,463. European Patent EP0593685. SafeSheath is a registered trademark of Pressure Products Medical Supplies, Inc. Made in USA