STERILE | EO

Pressure Products, Inc.
Customer Service
1861 N. Gaffey Street, Suite B
San Pedro, CA 90731 - USA
Tel +1-310-547-4973
Fax +1-310-547-4760
www.pressure-products.com

Oscor Inc.
3816 DeSoto Boulevard
Palm Harbor, FL 34683 - USA
Tel +1-727-937-2511
Fax +1-727-934-9835
www.oscor.com

EC | REP

Oscor Europe GmbH
Fritz-Vomfelde-Strasse 6
40547 Düsseldorf, Germany

CE
0344

SafeSheath is a registered trademark of Pressure Products, Inc.

**Pressure®
Products**
ACCESSING THE HEART OF INNOVATION

en Worley Right-sided Access Introducer System

cs Pravostraný přístupový zaváděcí systém Worley	1
da Worley Introducer-system med højre sides adgang	2
de Worley Einführungshilfe-System für Rechtskatherisator	3
el Σύστημα εισαγωγέα προστάλωσης & δεξιάς πλευράς Worley	3
es Sistema introducer con acceso lateral derecho Worley	4
et Eestikeelse Kasutusjuhend	5
fi Käytöohjeet	5
fr Consignes d'utilisation	6
hu Használati útmutató	6
it Istruzioni per l'uso	7
ja 使用説明書	7
ko 사용설명서	8
lt Naudojimo instrukcijos	9
nl Gebruiksaanwijzing	9
no Bruksanvisning	10
pl Instrukcja obsługi	10
pt Instruções de Utilização	11
ru Инструкции по использованию	11
sk Návod na použitie	12
sv Bruksanvisning	13
tr Kullanım Talimatı	13
zh 使用说明	14

en Instructions for Use	1
cs Návod k použití	2
da Brugsanvisning	2
de Gebrauchsanleitung	3
el Δείτε τις οδηγίες χρήσης	3
es Instrucciones de uso	4
et Eestikeelse Kasutusjuhend	5
fi Käytöohjeet	5
fr Consignes d'utilisation	6
hu Használati útmutató	6
it Istruzioni per l'uso	7
ja 使用説明書	7
ko 사용설명서	8
lt Naudojimo instrukcijos	9
nl Gebruiksaanwijzing	9
no Bruksanvisning	10
pl Instrukcja obsługi	10
pt Instruções de Utilização	11
ru Инструкции по использованию	11
sk Návod na použitie	12
sv Bruksanvisning	13
tr Kullanım Talimatı	13
zh 使用说明	14


Figure 1


Figure 2

en-Instructions for Use

This device is intended for one time use only. Read instructions prior to use.

Indications

For the introduction of various types of pacing leads and catheters.

Warnings

- This product is sensitive to light. Do not use if stored outside the protective outer carton. Store in a cool, dark, and dry place.
- Infusion through the side port can be done only after all air is removed from the unit.

Precautions

- Do not alter this device in any way.
- When assembling the sheath and dilator, care must be taken to insert the dilator tip straight through the center of the valve membrane in order to prevent inadvertent puncturing of the membrane.
- Aspiration and saline flushing of the sheath, dilator, and valve should be performed to help minimize the potential for air embolism and clot formation.
- Indwelling introducer sheaths should be internally supported by a catheter, electrode, or dilator.
- Dilators, catheters, and pacing leads should be removed slowly from the sheath. Rapid removal may damage the valve members resulting in blood flow through the valve. Never advance or withdraw guide wire or sheath when resistance is met. Determine cause by fluoroscopy and take remedial action.
- When injecting or aspirating through the sheath, use the side port only.
- **Single Use Devices:** This single-use product is not designed or validated to be reused. Reuse may cause a risk of cross-contamination, affect the measurement accuracy, system performance, or cause a malfunction as a result of the product being physically damaged due to cleaning, disinfection, re-sterilization, or reuse.

Use Sterile Technique A suggested procedure:

1. Peel open package and place contents on sterile field.
2. Prep skin and drape in area of anticipated puncture as desired.
3. The angle of the needle should be adjusted depending on the Patient's build: shallow in a thin person, deeper in a heavyset person.
4. Insert needle into vessel. The needle position should be verified by blood return.
5. Aspirate the puncture needle using the syringe.
6. Remove the syringe and insert soft tip of guide wire through the introducer needle into the vessel. Advance the guide wire to the required depth. Leave an appropriate amount of guide wire exposed.
7. Hold guide wire in place and remove introducer needle. Do not withdraw the guide wire back into the cannula as this may result in separation of the guide wire. The cannula should be removed first.
8. Insert vessel dilator into sheath and rotate the dilator cap over valve housing to secure the dilator onto sheath assembly.
9. Thread the dilator/sheath assembly over the guide wire.
10. Advance the dilator and sheath together with a twisting motion over the guide wire and into the vessel. Fluoroscopic observation may be advisable. Attaching a clamp or hemostat to the proximal end of the guide wire will prevent inadvertently advancing the guide wire entirely into the patient.
11. Once assembly is introduced into the vascular system, separate the dilator cap from the sheath valve housing by rotating the dilator cap off the hub. (see figure 1)
12. Slowly retract the guide wire and dilator, leaving the sheath in position. The hemostasis valve will reduce the loss of blood and the inadvertent aspiration of air through the sheath.
13. Aspirate all air from the sheath valve assembly by using a syringe connected to the side port. Flush the introducer through the side port. If the introducer is to remain in place during catheter positioning and testing, flushing the introducer via the side port periodically with saline is advised.
14. When the tip is positioned in the high right atrium, the hub of the sheath is designed to protrude 3-5 cm from the pocket depending on the size of the patient, where the sheath enters the vein, and the size of the RA.
15. For coronary sinus access, an appropriate sized access catheter may be introduced over the wire through the introducer valve, ie., a SafeSheath LVI Renal. See section titled **Coronary Sinus access using a SafeSheath LVI Renal with Target Vein Selector**.
16. Flush the sheath with 5cc of saline immediately before peeling the sheath away in order to minimize backbleeding.
17. First, withdraw the sheath from the vessel over the lead until the lead can be stabilized as it exits the muscle. Remove the sheath by sharply snapping the tabs of the valve housing in a plane perpendicular to the long axis of the sheath to split the valve and peel the sheath apart. (see figure 2)

Coronary Sinus access using a SafeSheath LVI Renal with Target Vein Selector (Model no. LVI/75-5-62-07-RE) A suggested procedure:

1. Position the Right-sided sheath with the tip in the high RA directed toward the coronary sinus. In this position, the proximal end (hub) of the sheath is typically 5-7cm out of the body.
2. Attach a contrast injection system to the side port of the SafeSheath® LVI Renal (LVI) and flush the LVI with contrast.
3. Advance the LVI through the hemostatic valve of the sheath.
4. Using fluoroscopic observation, advance the LVI until the tip of the LVI emerges from the tip of the sheath into the right atrium. A guide wire may be used through the LVI to aid this advancement.
5. Using fluoroscopic observation, apply clockwise torque and advance the LVI toward the RV.
6. Once the tip of the LVI has reached the RV or tricuspid annulus, apply counterclockwise torque to direct the tip of the LVI toward the coronary sinus ostium.

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