**For U.S.-California Only.**  
**Proposition 65**, a State of California voter initiative, requires the following notice:

**WARNING:** This product and its packaging have been sterilized with ethylene oxide. This packaging may expose you to ethylene oxide, a chemical known to the state of California to cause cancer or birth defects or other reproductive harm.

**Pressure Products, Inc.**  
Sales/Customer Service  
1861 N. Gaffey Street, Suite B  
San Pedro, CA 90731 - USA  
Tel +1-800-600-4973  
Fax +1-310-547-4760  
[www.pressure-products.com](http://www.pressure-products.com)

**Martech Medical Products**  
1500 Delp Drive  
Harleysville, PA 19438  
Tel +1-215-256-8833  
Tel +1-215-256-8837  
[www.martechmedical.com](http://www.martechmedical.com)

**EC REP**

**MPS Medical Product Service GmbH**  
Bornsgasse 20  
35619 Braunschweig, Germany  
Tel +49-(0)-6442-962073  
Fax +49-(0)-6442-962074  
[info@mps-gmbh.eu](mailto:info@mps-gmbh.eu)


**CE**

**0434**

**CAUTION: Federal (U.S.A.) Law restricts this device to sale by or on the order of a physician.**

## Technical Manual

<b>en</b>	<b>Directions for Use</b>	<b>1</b>
<b>cs</b>	<b>Návod k použití</b>	<b>2</b>
<b>da</b>	<b>Brugsvejledning</b>	<b>2</b>
<b>de</b>	<b>Gebrauchsanweisung</b>	<b>3</b>
<b>el</b>	<b>Οδηγίες χρήσης</b>	<b>3</b>
<b>es</b>	<b>Instrucciones de uso</b>	<b>4</b>
<b>et</b>	<b>Kasutusjuhend</b>	<b>4</b>
<b>fi</b>	<b>Käyttöohje</b>	<b>5</b>
<b>fr</b>	<b>Mode d'emploi</b>	<b>5</b>
<b>hu</b>	<b>Használati útmutató</b>	<b>6</b>
<b>it</b>	<b>Istruzioni per l'uso</b>	<b>6</b>
<b>ja</b>	<b>使用説明書</b>	<b>7</b>
<b>ko</b>	<b>사용 지침</b>	<b>7</b>
<b>lt</b>	<b>Naudojimo instrukcijos</b>	<b>8</b>
<b>nl</b>	<b>Gebruiksaanwijzing</b>	<b>8</b>
<b>no</b>	<b>Bruksanvisning</b>	<b>9</b>
<b>pl</b>	<b>Instrukcja użycia</b>	<b>9</b>
<b>pt</b>	<b>Instruções de Utilização</b>	<b>10</b>
<b>ru</b>	<b>Инструкция по применению</b>	<b>11</b>
<b>sk</b>	<b>Pokyn na použitie</b>	<b>11</b>
<b>sv</b>	<b>Bruksanvisning</b>	<b>12</b>
<b>tr</b>	<b>Kullanım Talimatları</b>	<b>12</b>
<b>zh</b>	<b>使用说明</b>	<b>13</b>

## en-Directions for Use

The SafeVue™ Balloon Catheter is a dual lumen catheter for temporary vascular occlusion or balloon catheter flotation while providing a port for distal contrast injection or guidewire placement. The catheter may be used for procedures requiring angiography, wedge pressure measurements, or any other procedure where a vessel needs to be occluded while maintaining a distal lumen.

### Indications

For temporary vascular occlusion and balloon catheter flotation with contrast infusion or guidewire placement.

### Warnings

- For one use only, do not re-sterilize.
- Store in a dry, dark, cool place.
- Do not use if package is open or damaged.
- Filtered CO<sub>2</sub> must be used to inflate the balloon if there is a possibility that balloon rupture would result in air embolus in the heart or arterial circulation. Liquids must not be used as a balloon inflation medium.
- Although diagnostic cardiac catheterization procedures have been proven to be safe, the physician should be aware that certain complications can occur and it is recommended that the user of this product become familiar with the guidelines established by Drs. Swan and Ganz for the safe use of the SafeVue™ Balloon Catheter. (Swan HJC, Ganz W, Forrester, et al: Catheterization of the heart in man with the use of a flow directed balloon tipped catheter. *N Engl J Med* 283:447-451, 1970.)

### Precautions

- Precautions
- Do not alter this device in any way.
- Aspiration and flushing of the lumen should be performed to help minimize the potential for air embolism and clot formation.
- A continuous drip or an intermittent flush should be used to maintain patency of the lumen.
- When using the Balloon Catheter with a valved introducer, place the Transvalvular Insertion tool (TVI) through the introducer hemostasis valve prior to insertion or withdrawal of the Balloon Catheter.
- Always deflate the balloon by pulling back on the syringe prior to inserting or withdrawing the Balloon Catheter.
- Do not advance or retract Balloon Catheter if resistance is met. Determine the cause of resistance before proceeding.

### Catheter Inspection and Testing

The SafeVue™ Balloon Catheter incorporates a pressure relief valve in fluid communication with the balloon inflation lumen to minimize the potential for vessel trauma or dissection. Inflation of the balloon higher than the set pressure of the relief valve will divert the excess air (or CO<sub>2</sub> inflation medium) to escape from the hub.

1. Remove the Balloon Catheter from the package using aseptic technique.
2. Slide the Transvalvular Insertion tool (TVI) towards the balloon catheter hub to completely expose the balloon.
3. Inflate the balloon to capacity and immerse the balloon in sterile room temperature saline. If there is any evidence of air bubbles escaping the balloon or if the balloon will not remain inflated, do not use the Balloon Catheter.
4. At no time should liquid be used for balloon inflation.
5. Aspirate the air completely from the balloon and close the stopcock to prevent the balloon from opening.
6. Flush the inner catheter lumen with sterile heparinized saline.
7. To remove the Balloon Catheter, open the stopcock, and with a syringe connected to the stopcock, pull back on the syringe until resistance is met. Close the stopcock and retract the Balloon Catheter. When using a valved introducer, place the Transvalvular Insertion tool (TVI) through the hemostasis valve prior to withdrawing the Balloon Catheter.

### Use Sterile Technique A suggested procedure:

Techniques may vary according to physician preference or the patient's anatomy or physical condition.

1. After inspecting and testing the SafeVue™ Balloon Catheter, it may be introduced into the body by inserting it through an appropriately sized introducer or guide catheter. When using a valved introducer, place the Transvalvular Insertion tool (TVI) through the hemostasis valve prior to inserting the Balloon Catheter (*Figure 1*).


**Always deflate the balloon by pulling back on the syringe prior to inserting or withdrawing the balloon catheter. Do not advance or retract Balloon Catheter if resistance is met. Determine the cause of resistance before proceeding.**

Fluoroscopic verification of the Balloon Catheter's position is suggested.

2. Depth markings are provided along the catheter shaft at 10cm increments.
3. Once the Balloon Catheter is in the desired position, fill the syringe with 3.0 cc of inflation medium (air or filtered CO<sub>2</sub>).
4. Attached the syringe to the inflation lumen, open the stopcock, and slowly inject the inflation medium into the inflation lumen to inflate the balloon. Stop inflating the balloon when resistance is felt or a vibrating sound is heard at the hub indicating that the balloon relief pressure has been met. Close the stopcock.
5. Aspirate blood through the infusion lumen to remove air from the lumen.
6. With the inflated balloon occluding the vessel, contrast may be injected through the infusion lumen or a guidewire may be placed.

7. To remove the Balloon Catheter, open the stopcock, and with a syringe connected to the stopcock, pull back on the syringe until resistance is met. Close the stopcock and retract the Balloon Catheter. When using a valved introducer, place the Transvalvular Insertion tool (TVI) through the hemostasis valve prior to withdrawing the Balloon Catheter.

**Do not advance or retract Balloon Catheter if resistance is met. Determine the cause of resistance before proceeding.**

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