
SafeSept®


Transseptal Guidewire

TECHNICAL MANUAL

en	Instructions for Use	2
ar	تعليمات الاستخدام	2
bg	Инструкции за употреба	3
cs	Návod k použití	3
da	Brugsanvisning	4
de	Gebrauchsanleitung	4
el	Δείτε τις οδηγίες χρήσης	5
es	Instrucciones de uso	6
et	Eestikeelne Kasutusjuhend	6
fi	Käyttöohjeet	7
fr	Instructions d'utilisation	8
hr	Upute za uporabu	8
hu	Használati útmutató	9
it	Istruzioni per l'uso	9
ja	使用説明書	10
kk	Пайдалану жөніндегі нұсқаулық	10
ko	사용설명서	11
lt	Naudojimo instrukcijos	12
lv	Lietošanas pamācība	12
nl	Gebruiksaanwijzing	13
no	Bruksanvisning	13
pl	Instrukcja użycia	14
pt (br)	Instruções de uso	15
pt (ib)	Instruções de Utilização	15
ro	Instrucțiuni de utilizare	16
ru	Инструкция по применению	16
sl	Navodila za uporabo	17
sk	Návod na použitie	17
sr	Uputstvo za upotrebu	18
sv	Bruksanvisning	19
tr	Kullanma talimatı	19
zh	使用说明	20

EC REP

Bisping Medizintechnik GmbH
An der Glashütte 8
52074 Aachen - Germany
Tel +49 (0) 241-173518
Fax +49 (0) 241-175627

Pressure Products, Inc.
Customer Service
1861 N. Gaffey Street, Suite B
San Pedro, CA 90731 - USA
Tel +1-310-547-4973
Fax +1-310-547-4760
www.pressure-products.com

PPMDM
1 School Street
Morton, PA 19070 - USA
Tel +1-610-285-9858
Fax +1-610-285-9859

Rx ONLY CAUTION: Federal (U.S.A.) Law restricts this device to sale by or on the order of a physician.

SafeSept is a registered trademark of Pressure Products Medical Supplies, Inc. Patents 9,821,145; 9,585,692; 8,992,556; 8,500,697; 8,292,910; 8,157,829; 7,666,203; EP1542593; EP2259832; EP2647405; EP2668972; JP4496223. Other U.S. and worldwide patents pending. Protected by Patent Insurance.

en-Instructions for Use

The SafeSept Transseptal Guidewire® is a nitinol guidewire designed to cross the interatrial septum when supported by a transseptal introducer and be atraumatic when advanced unsupported into the left atrium. A radiopaque coil is positioned on the distal end of the device to provide fluoroscopic localization during the procedure. White printed markers are positioned on the proximal end of the device to provide approximate Transseptal Guidewire tip location relative to the needle tip during the procedure for use either with or without a hemostatic Y-adaptor connected to the needle hub.

INTENDED USE

The SafeSept Transseptal Guidewire® is used in conjunction with a transseptal needle to create the primary puncture in the interatrial septum and to guide the needle, dilator, and introducer through the septum from the right side of the heart to the left side.

INDICATIONS FOR USE

SafeSept® is indicated for use in procedures where access to the left atrium via the transseptal technique is desired. The SafeSept Transseptal Guidewire® is intended for single use only.

CONTRAINDICATIONS

The use of the Transseptal Guidewire is contraindicated in patients with the following conditions:

- Distorted anatomy due to congenital heart disease or other causes
- Significant chest or spine deformity
- The inability to lie flat
- Ongoing anticoagulation
- Left atrial thrombus or tumor
- Dilated aortic root
- Previous patch repair of the interatrial septum

WARNINGS AND PRECAUTIONS

- Single-Use Only: Do not re-use this device. After use thorough cleaning of biological and foreign material is not possible. Adverse patient reactions may result from re-use of this device.
- Store in a cool, dark, and dry place
- Prolonged exposure to temperatures above 30°C (86 °F) may damage the product
- Do not use if package is open or damaged in any way
- Do not use if tip protector not in place

ADVERSE EVENTS

In addition to all of the complications associated with any transseptal cardiac catheterization, the following can occur during the use of the Transseptal Guidewire:

- Puncture of the atrial free wall
- Puncture of the aorta
- Puncture of the inferior vena cava
- Puncture of the coronary sinus
- Tamponade
- Hemothorax
- Arterial embolism from thrombus at the puncture site
- Pulmonary embolism
- Stroke
- Death
- Atrial arrhythmias
- Residual atrial septal defects

INSTRUCTIONS FOR USE

Follow standard transseptal technique:

1. Per standard technique, aspirate and flush the transseptal system (including the needle, dilator, and sheath) before and after it is inserted into the right atrium.
2. Remove the protective tip cover and back the Transseptal Guidewire slowly into the guidewire tip straightener. Do not pull the Transseptal Guidewire completely out of the tip straightener.
3. Fully insert the tip straightener into the hub of the transseptal needle.
4. Advance the Transseptal Guidewire until the tip enters the curved portion of the transseptal needle (Printed markers are provided as a guide for common commercially available needle systems. The Transseptal Guidewire can generally be advanced to the first printed marker when using a transseptal needle without a Y-adaptor. Do not advance the guidewire beyond the tip of the needle).
5. Tent the fossa ovalis with the transseptal dilator. The transseptal needle tip position should be less than 5cm back from the tip of the dilator. يكون وراء رأس الإبرة العابرة للحاجز أقل من 0,5 سم وراء رأس الموسع.
6. Under visual guidance, while tenting and maintaining constant force on the septum with the dilator, slowly advance the Transseptal Guidewire through the transseptal needle, dilator, and across the interatrial septum, into the left atrium. Continue to advance into one of the pulmonary veins (The radio-opaque coil on the guidewire should be seen to be within the left atrium and subsequently the pulmonary vein). **At no time should the guide wire be advanced or withdrawn when resistance is met without first determining the cause visually and taking remedial action.**
7. The Transseptal Guidewire should be advanced approximately 7cm into the left atrium so that the entire radiopaque coil is across the septum and visible to ensure that the transseptal needle tip does not engage the radiopaque coil during advancement.
8. Slowly advance the transseptal needle through the dilator and over the Transseptal Guidewire across the septum. **Do not hold or pin the Transseptal Guidewire when advancing the transseptal needle.** If desired, measure pressures and/or inject contrast through the needle before advancing the dilator and introducer sheath. Once the needle tip position is confirmed in the left atrium, advance the dilator and introducer sheath over the Transseptal Guidewire into the left atrium.
9. Slowly remove the Transseptal Guidewire, needle, and dilator as a unit, leaving the introducer sheath in place. **Do not remove the Transseptal Guidewire first, leaving the transseptal needle point exposed inside the left atrium. At no time should the introducer, needle, or guide wire be advanced or withdrawn when resistance is met without first determining the cause visually and taking remedial action.**
10. After use, the Transseptal Guidewire may be a potential biohazard. Handle and dispose of it in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

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