# Transvalvular Insertion Tool (TVI)

STERILE EO

**CAUTION: Federal (U.S.A.) Law restricts this device to sale by or on the order of a physician.**

**Pressure Products**  
Customer Service  
1861 N. Gaffey Street, Suite B  
San Pedro, CA 90731 - USA  
Tel +1-310-547-4973  
Fax +1-310-547-4760  
www.pressure-products.com

**Merit Medical Systems, Inc.**  
65 Great Valley Parkway  
Malvern, PA 19355 - USA  
Tel +1-801-253-1600  
USA Customer Service:  
Tel +1-800-356-3748


**Merit Medical Ireland, Ltd**  
Parkmore Business Park West  
Galway, Ireland  
European Customer Service:  
Tel +31 43 358 82 22

0344

## Technical Manual

en	Directions for Use	1
cs	Návod k použití	1
da	Brugsanvisning	1
de	Gebrauchsanleitung	2
el	Οδηγίες χρήσης	2
es	Instrucciones para el uso	2
et	Kasutusjuhend	2
fi	Käyttöohjeet	3
fr	Mode d'emploi	3
hu	Használati utasítás	3
it	Istruzioni per l'uso	3
ja	使用説明書	4
ko	사용방법	4
lt	Naudojimo nurodymai	4
nl	Gebruiksaanwijzing	4
no	Bruksanvisning	5
pl	Instrukcja użytkowania	5
pt	Instruções de Utilização	5
ru	Инструкция по применению	5
sk	Pokyny na používanie	6
sv	Bruksanvisning	6
tr	Kullanım Talimatları	6
zh	使用方法	6

## en-Directions for Use

This device is intended for one time use only. Read instructions prior to use.

### Indications

For the introduction of various types of pacing or defibrillator leads and catheters.

### Warnings

Improper use of the transvalvular insertion tool (TVI) can cause air embolism and back bleeding.

### Precautions

- Do not alter this device in any way.
- Federal law (U.S.A.) restricts this device to sale by or on the order of a physician.
- When using the transvalvular insertion tool (TVI), lead size may not exceed 6.2F for the 7F TVI, or 9.0F for the 9F TVI.
- When using the TVI always keep the exposed proximal end covered to prevent air embolization and back bleeding.
- **Single Use Devices:** This single-use product is not designed or validated to be reused. Reuse may cause a risk of cross-contamination, affect the measurement accuracy, system performance, or cause a malfunction as a result of the product being physically damaged due to cleaning, disinfection, re-sterilization, or reuse.

### Use Sterile Technique *A suggested procedure:*

1. A transvalvular insertion tool (TVI) is to be used at the physician's discretion to open the valve for ease of lead placement.
2. **CAUTION! When the TVI is inserted into the SafeSheath valve housing all hemostasis is lost and the risk of air embolization and back bleeding exists. Always keep the proximal exposed end of the TVI covered with your thumb while the TVI is in use.**
3. To use the TVI, insert the distal end of the TVI into the valve housing by gently pushing the TVI into the sheath.
4. Hold your thumb over the proximal exposed opening of the TVI to prevent air embolization or back bleeding.
5. Advance the pacemaker lead through the TVI and into the sheath.
6. As soon as the pacer lead is resting inside the sheath pull the TVI back out of the sheath valve housing.
7. The TVI may then either be peeled away or it may be temporarily left resting on the shaft of the pacemaker lead.
8. Once the TVI has been withdrawn from the sheath, aspirate the sheath, through the sideport, until any air which may have entered the sheath during the procedure is removed and again flush with Heparinized saline.

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