Transvalvular Insertion Tool (TVI)

STERILE EO

CAUTION: Federal (U.S.A.) Law restricts this device to sale by or on the order of a physician.

Pressure Products, Inc.
Customer Service
1861 N. Gaffey Street, Suite B
San Pedro, CA 90731 - USA
Tel +1-800-600-4973
Tel +1-310-547-4760
www.pressure-products.com

Thomas Medical Products
A GE Healthcare Company
65 Great Valley Parkway
Malvern, PA 19355 - USA
Tel +1-610-651-5000
Tel +1-610-651-5003
Fax +1-610-651-5003

GE Medical Systems SCS
283 rue de la Minière
78530 BUC, FRANCE
Tel +33 130704040
Fax +33 130704440

0044

Technical Manual

en-Directions for Use

This device is intended for one time use only. Read instructions prior to use.

Indications

For the introduction of various types of pacing or defibrillator leads and catheters.

Warnings

Improper use of the transvalvular insertion tool (TVI) can cause air embolism and back bleeding.

Precautions

- Do not alter this device in any way.
- Federal law (U.S.A.) restricts this device to sale by or on the order of a physician.
- When using the transvalvular insertion tool (TVI), lead size may not exceed 6.2F for the 7F TVI, or 9.0F for the 9F TVI.
- When using the TVI always keep the exposed proximal end covered to prevent air embolization and back bleeding.
- Single Use Devices:** This single-use product is not designed or validated to be reused. Reuse may cause a risk of cross-contamination, affect the measurement accuracy, system performance, or cause a malfunction as a result of the product being physically damaged due to cleaning, disinfection, re-sterilization, or reuse.

Use Sterile Technique *A suggested procedure:*

- A transvalvular insertion tool (TVI) is to be used at the physician's discretion to open the valve for ease of lead placement.
- CAUTION! When the TVI is inserted into the SafeSheath valve housing all hemostasis is lost and the risk of air embolization and back bleeding exists. Always keep the proximal exposed end of the TVI covered with your thumb while the TVI is in use.**
- To use the TVI, insert the distal end of the TVI into the valve housing by gently pushing the TVI into the sheath.
- Hold your thumb over the proximal exposed opening of the TVI to prevent air embolization or back bleeding.
- Advance the pacemaker lead through the TVI and into the sheath.
- As soon as the pacer lead is resting inside the sheath pull the TVI back out of the sheath valve housing.
- The TVI may then either be peeled away or it may be temporarily left resting on the shaft of the pacemaker lead.
- Once the TVI has been withdrawn from the sheath, aspirate the sheath, through the sideport, until any air which may have entered the sheath during the procedure is removed and again flush with Heparinized saline.

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