

Coronary Sinus **Balloon Venography Catheter**

An aid for left ventricular epicardial lead placement.

The Pressure Products Coronary Sinus **Balloon Venography Catheter** is designed for placement in the proximal coronary sinus via the SafeSheath CSG introducer. When inflated, the coronary sinus is occluded and angiography easily accomplished.

The Balloon Venography Catheter is part of a complete line of implantable device therapy products and accessories available from Pressure Products.

Each package contains:

- 1cm diameter Balloon
- 2 Lumens
- Soft J Tip—6 French, 80cm long
- Volume-limited Syringe

Model	Description	Size	Length
BVCS 6180	Coronary Sinus Balloon Venography Catheter	6F	80cm

Pressure Products, Inc.

Customer Service

1861 N. Gaffey Street Suite B San Pedro, CA 90731

(800) 600-4973 (310) 547-4973 Fax (310) 547-4760 www.pressure-products.com