


SiteLIGHT™

STERILE SURGICAL LIGHT

Provides a cool, bright, focused white light exactly where needed during surgical procedures.



■ Super bright white LED with 5100° Kelvin color temperature

■ Top-mounted pushbutton on/off switch

Introducing the SiteLight™ disposable sterile surgical light. The super small design of the SiteLight enables you to put a focused beam of cool white light directly where you need it during surgical procedures.

Incorporating a convenient top-mounted pushbutton switch, the light can be turned on and off as needed, providing over 60 minutes of continuous operation. And with a color temperature of approximately 5100° K, the SiteLight provides virtually tint-free illumination of the site for accurate viewing.


The self-adhesive disk enables the SiteLight to be releasably mounted to any sterile surface in the surgical setting, including gloves, drapes, surgical tools or other sterile devices.

The SiteLight disposable sterile surgical light is part of a complete

line of implantable device therapy products and accessories from Pressure Products.

Features:

- Sterile, low-cost disposable design
- Super bright white light with approximately 5100° Kelvin color temperature
- Very light, weighing less than 9 grams
- Compact, low-profile design with self-adhesive pad for easy mounting to any surface.
- Fixed-focus 2 inch diameter LED beam
- Over 60 minutes of continuous light output
- Cool to the touch



■ Self-adhesive, easily releasable pad allows mounting to any convenient surface or device

Model No.	Description	Color Temperature
SSL	Disposable sterile surgical light	5100° Kelvin LED

USA and Worldwide Patents pending.