


# Transvalvular Insertion Tool (TVI)

*An aid to lead placement when using the SafeSheath® family of introducers*




The Pressure Products **Transvalvular Insertion Tool (TVI)** is designed to open the valve for lead placement via the SafeSheath family of introducers. After lead placement, the TVI is easily peeled away.

The Transvalvular Insertion Tool is part of a complete line of implantable device therapy products and accessories available from Pressure Products.

### Features:

- Aids in opening valve for ease of lead placement
- Can be used with leads up to 9 French
- Is easily peeled away after lead placement



Model	Description	Largest Lead
TVI-07	Transvalvular Insertion Tool (TVI)	6.2F
TVI-09	Transvalvular Insertion Tool (TVI)	9F