


# Transvalvular Insertion Tool (TVI)




*An aid to lead placement when using the SafeSheath® family of introducers*

The Pressure Products **Transvalvular Insertion Tool (TVI)** is designed to open the valve for lead placement via the SafeSheath family of introducers. After lead placement, the TVI is easily peeled away.

The Transvalvular Insertion Tool is part of a complete line of implantable device therapy products and accessories available from Pressure Products.

### **Features:**

- Aids in opening valve for ease of lead placement
- Can be used with leads up to 6.2 French
- Is easily peeled away after lead placement



| Model  | Description                        | Largest Lead |
|--------|------------------------------------|--------------|
| TVI-07 | Transvalvular Insertion Tool (TVI) | 6.2F         |